

Specification Sheet



Description: Corrosive Cabinet - 60L

Product Code: 04-1071

Corrosive Cabinet - 60L

Product Features

A taller version of the 30L cabinet, this model is a critical storage unit for smaller containers in laboratories, workshops and areas with limited corrosive requirements.

- Built to comply with the requirements of AS3780
- Lockable, non-sparking 3 point latching
- Self-closing door system
- Double walled 1.2mm zinc steel construction
- Fitted with ventilation ducts and built in flame arrestors
- Includes Class 8 hazard labelling and other decals required by the standard
- 150mm liquid tight sump
- Sturdy galvanized base shelf that can be removed for easy cleaning
- Adjustable shelf slots – comes with one shelf as standard (additional shelves available)



chemshed

External Dimensions (w,d,h)	Internal Dimensions (w,d,h)	Capacity	No. of Shelves	Weight
500 x 450 x 1,065mm	418 x 370 x 818mm	60L	1 (plus base)	62kgs

To request a quote or need further information, contact us now!
Call 0800 951 010 or email customerservicenz@dksh.com