

Specification Sheet



Description: Corrosive Cabinet - 100L

Product Code: 04-1072

Corrosive Cabinet - 100L

Product Features

A squat, double door cabinet that is designed to fit under a bench and securely store a variety of size containers of Class 8 corrosive liquids.

- Built to comply with the requirements of AS3780
- Lockable, non-sparking 3 point latching
- Self-closing door system
- Double walled 1.2mm zinc steel construction
- Fitted with ventilation ducts and built in flame arrestors
- Includes Class 8 hazard labelling and other decals required by the standard
- 150mm liquid tight sump
- Sturdy galvanized base shelf that can be removed for easy cleaning



chemshed

External Dimensions (w,d,h)	Internal Dimensions (w,d,h)	Capacity	No. of Shelves	Weight
920 x 615 x 800mm	838 x 535 x 553mm	100L	1 (plus base)	100kgs

**To request a quote or need further information, contact us now!
Call 0800 951 010 or email customerservicenz@dksh.com**