

# Specification Sheet



**Description:** Toxic Cabinet - 30L

**Product Code:** 04-1140

## Toxic Cabinet - 30L

### Product Features

A Chemshed under bench storage cabinet for toxic substances with a 30L capacity and featuring the patented self closing system.

A solid and reliable cabinet for storing smaller quantities of Class 6 Toxic substances in a safe and compliant manner. Comes with the required Class 6 hazard labelling, this cabinet is built to exceed the latest design requirements for complaint storage.

- Built to comply with the requirements of AS/NZ4452 standard
- Self-closing, non-latching door (as required)
- Built in flash arrestors and vent openings
- High capacity, liquid-tight sump
- Supplied with Class 6 decals



chemshed

External Dimensions (w,d,h)	Internal Dimensions (w,d,h)	Capacity	No. of Shelves	Weight
500 x 450 x 800mm	418 x 370 x 553mm	30L	1 (plus base)	33kgs

**To request a quote or need further information, contact us now!**  
**Call 0800 951 010 or email [customerservicenz@dksh.com](mailto:customerservicenz@dksh.com)**