

## Universal ovens

### AP series large capacity

#### Applications

- Drying.
- Set temperature.
- Stress tests.
- Sterilization.
- The series of ovens AP was designed for all thermal testing precision room temperature + 10 ° C to + 250 ° C.

#### Spécifications

- 2 volumes 414 and 714 liters.
- Stainless steel interior.
- Models equipped with two perforated shelves in anti-tip safety.
- Display simultaneously at 1/10 ° C setpoint and actual temperature.
- Uniformity at 150 ° C +/- 1.7 ° C
- Control by Pt 100 Ohms.
- Double Door leaves.
- Mounted on castors.
- IP protection to EN 60529.
- Rated power: 3100 W.



Reference	Description	Ext. dimensions W x D x H (mm)	Int. Dimensions W x D x H (mm)	Nb of shelves
AP414/230	Universal Oven Model FI407U 407 liters 380 V tri + N	1120 x 700 x 1160	980 x 500 x 740	2
AP414	Universal Oven Model FI407U 407 liters 230 V tri	1120 x 700 x 1160	980 x 500 x 740	2
AP7143	Universal Oven Model FI407U 407 liters 380 V tri	1120 x 700 x 1160	980 x 500 x 740	2
AP714	Universal Oven Model FI727U 727 liters 380 V tri + N	1140 x 800 x 1540	1000 x 600 x 1140	2

#### ACCESSORIES

Reference	Description
AP/ETA3	Additional stainless steel shelf
FI/OPT1	Air extraction adjustable flow from 0 to 200m3/hour
AP/OPT7	
FI/OPT4	
AP/OPT8	
FI/OPT8	

#### ABOUT US

Founded in Lyon in 1918 at the request of pharmaceutical companies which then grow in the Lyon region and looking for the sample preparation equipment...

#### SERVICE DEPARTMENT

From 8.00 to 17.00 Monday to Thursday  
From 8.00 to 16.00 on Fridays

#### CONTACT OUR SERVICE

+33 (0)1 60 95 15 70

## CONTACT US

Froilabo

5 avenue Lionel Terray 69330 Meyzieu - FRANCE

Pho : +33 (0)4 78 04 75 75 - Fax : +33 (0)4 78 93 08  
24

35 boulevard Beaubourg 77184 Emerainville -  
FRANCE

Pho : +33 (0)1 60 05 15 65 - Fax : +33 (0)1 60 77 41

[Sitemap](#) - [Legal notices](#)