



Incubators and refrigerated incubators

Static series BSP

Applications

- Natural convection is designed for long duration incubation, incubating substrates sensitive to dryness and also powder applications.



Specifications

- Combined PID / fuzzy logic regulation
- Pt 100 Ohms control sensor
- Digital Display 1/10° accurate (actual and set-point)
- Timer: 99h 59 min
- High and low temperature visual alarm
- Over-temperature safety according to the European standard, NF EN61010-2-010.
- Pre-heated chamber (accurate temperature control preventing sample desiccation)
- Conforms to standards, EN60068-3-11 and NFX15-140
- Easy "hands-free" closing door
- EASYCLEANING inner stainless steel chamber (flat surface, removable shelf holders)
- Individual quality control certificate
- 24 month warranty



Reference	Description	Voltage	Ext. dimensions W x D x H (mm)	Int. Dimensions W x D x H (mm)	Nb of shelves	Shelf Dimensions W x D (mm)
BSP65	Natural convection incubator, plain door 61 liters	220V, 50Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP65PV	Natural convection incubator, glass door 61 liters	220V, 50Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP65/26	Natural convection incubator, plain door 61 liters	220V, 60Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP65PV/26	Natural convection incubator, glass door 61 liters	220V, 60Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP65/156	Natural convection incubator, plain door 61 liters	110V, 50/60Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP65PV/156	Natural convection incubator, glass door 61 liters	110V, 50/60Hz	526 x 579.5 x 640	400 x 394 x 390	1	380 x 320
BSP125	Natural convection incubator, plain door 124 liters	220V, 50Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP125PV	Natural convection incubator, glass door 124 liters	220V, 50Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP125/26	Natural convection incubator, plain door 124 liters	220V, 60Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP125PV/26	Natural convection incubator, glass door 124 liters	220V, 60Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP125/156	Natural convection incubator, plain door 124 liters	110V, 50/60Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP125PV/156	Natural convection incubator, glass door 124 liters	110V, 50/60Hz	626 x 679.5 x 750	500 x 494 x 500	2	480 x 430
BSP245	Natural convection incubator, plain door 242 liters	220V, 50Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430
BSP245PV	Natural convection incubator, glass door 242 liters	220V, 50Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430
BSP245/26	Natural convection incubator, plain door 242 liters	220V, 60Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430
BSP245PV/26	Natural convection incubator, glass door 242 liters	220V, 60Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430
BSP245/156	Natural convection incubator, plain door 242 liters	110V, 50/60Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430
BSP245PV/156	Natural convection incubator, glass door 242 liters	110V, 50/60Hz	626 x 679.5 x 1230	500 x 494 x 980	2	480 x 430

ACCESSORIES

Reference	Model	Description
AP/RAY1	BSP65	Stainless steel shelf
AP/RAY2	BSP125 BSP245	Stainless steel shelf
BP/ETABP1	BSP65	Stainless steel wire rack for 6 x 6 ? 90 mm Petri boxes
BP/ETABP2	BSP125 BSP245	Stainless steel wire rack for 9 x 6 ? 90 mm Petri boxes
AP/SUP	All	Stacking kit

OPTIONS

Reference	Description
BP/OPT1	Cable entry port 50 mm
AP/OPT3	Key locking door
AP/OPT6	Fitted with 4 castors
AP/OPT8	Certified Pt 100 Ohms + DIN socket for recorder link
AP/OPT10	Circular 1 way disc recorder
AP/OPT11	Electronic recorder 1 way continuous graph
QUAL4	9 points temperature monitoring
QUAL3	Certificate for 1 oven - 1 set-point

ABOUT US

Founded in Lyon in 1918 at the request of pharmaceutical companies which then grow in the Lyon region and looking for the sample preparation equipment...

SERVICE DEPARTMENT

From 8.00 to 17.00 Monday to Thursday
From 8.00 to 16.00 on Fridays

CONTACT OUR SERVICE

+33 (0)1 60 95 15 70

CONTACT US

Froilabo

5 avenue Lionel Terray 69330 Meyzieu - FRANCE

Pho : +33 (0)4 78 04 75 75 - Fax : +33 (0)4 78 93 08
24

35 boulevard Beaubourg 77184 Emerainville -
FRANCE

Pho : +33 (0)1 60 95 15 70 - Fax : +33 (0)1 60 77 44

[Sitemap](#) - [Legal notices](#)