



## Incubators and refrigerated incubators

### BP Series high volume

#### Applications

- The FIB incubator range was developed to meet the requirements of biological and medical analysis laboratories, allowing them to have surpass the present standard.

#### Spécifications

- Temperature range: from ambient +7°C to +65°C.
- Ventilation to optimize temperature homogeneity
- Accurate Pt 100 Ohms probe controls
- Stainless steel interior
- Power supply 230V / 50Hz



Reference	Volume (liters)	Ext. Dimensions W x D x H (mm)	Int. Dimensions W x D x H (mm)
BP414	407	1120 x 700 x 1160	980 x 500 x 740
BP714	727	1140 x 800 x 1540	1000 x 600 x 1140

#### ACCESSORIES

Reference	Description
AP/ETA3	Stainless steel shelf

#### OPTIONS

Reference	Description
FI/OPT1	Puller has variable flow
AP/OPT7	50mm cable entry
FI/OPT4	Programmable controller 16 segments
AP/OPT8	Type of sensor PT100 ohm + socket din
FI/OPT8	Overtemperature audible alarm

#### ABOUT US

Founded in Lyon in 1918 at the request of pharmaceutical companies which then grow in the Lyon region and looking for the sample preparation equipment...

#### SERVICE DEPARTMENT

From 8.00 to 17.00 Monday to Thursday  
From 8.00 to 16.00 on Fridays

#### CONTACT OUR SERVICE

+33 (0)1 60 95 15 70

#### CONTACT US

Froilabo  
5 avenue Lionel Terray 69330 Meyzieu - FRANCE  
Pho : +33 (0)4 78 04 75 75 - Fax : +33 (0)4 78 93 08 24  
35 boulevard Beaubourg 77184 Emerainville - FRANCE  
Pho : +33 (0)1 60 95 15 70 - Fax : +33 (0)1 60 77 41 41

