



## Incubators and refrigertaed incubators

# Refrigerated high volume series - COOLSTORE

### Applications

- The "Coolstore" incubators are designed for germination or growth applications.
- The two temperature range options of 0°C and -20°C allow for more applications without affecting homogeneity or stability of the metrological performances.

### Spécifications

- Temperature range: 0°C to +60°C or -20°C to +60°C.
- Ventilation to optimize temperature homogeneity
- Stainless steel interior for easy cleaning and disinfection
- Key locking door
- Mounted on 4 castors, 2 with breaks
- Refrigerant without CFC/HCFC: R134a

### Avantages

- Efficient electronic PID controls
- Dual digital display 1/10° real and set-point
- Pt 100 Ohms probe
- Class 3.1 safety thermostat for over-temperature
- Cable entry port diam. 35 mm with lid



Volume (liters)	Voltage	Power	Ext. dimensions W x D x H (mm)	Int. Dimensions W x D x H (mm)	Weight (kg)	Nb of Shelves
600	220V/50Hz	1650 W	700 x 830 x 2130	530 x 660 x 1550	170 kg	4
1200	220V/50Hz	2400 W	1430 x 830 x 2130	1270 x 660 x 1550	210 Kg	2x4

### ACCESSORIES

Reference	Model	Description
MEDI600/ETA	Coolstore 600	Additional stainless steel shelf 1 shelf per level, W x D : 530x655 mm (max. weight 30 kg)
MEDI1300/ETA	Coolstore 1300	Additional stainless steel shelf 3 shelves per level, W x D : 530x655 mm (max. weight 30 kg)

### OPTIONS

Reference	Models	Description
BP/OPT1	All	Cable entry port ? 50 mm
AP/OPT8	All	Certified Pt 100 Ohms + DIN socket for recorder link
AP/OPT10	All	Circular 1 way disc recorder
AP/OPT11	All	Electronic recorder 1 way continuous graph
FI/OPT8	All	Sound alarm for over-temperature
COOLS/OPT1	All	Automatic defrost (recommended between 0 and +10°C)
COOLS/LIGHT1	Coolstore 600	Daily programmable shelf lighting 3 lighting levels
COOLS/LIGHT2	Coolstore 1300	Daily programmable shelf lighting 2 x 3 lighting levels
QUAL4	All	9 points temperature monitoring
QUAL3	All	Certificate for 1 oven - 1 set-point

## ABOUT US

Founded in Lyon in 1918 at the request of pharmaceutical companies which then grow in the Lyon region and looking for the sample preparation equipment...

## SERVICE DEPARTMENT

From 8.00 to 17.00 Monday to Thursday  
From 8.00 to 16.00 on Fridays

## CONTACT OUR SERVICE

+33 (0)1 60 95 15 70

## CONTACT US

Froilabo

5 avenue Lionel Terray 69330 Meyzieu - FRANCE

Pho : +33 (0)4 78 04 75 75 - Fax : +33 (0)4 78 93 08

24

35 boulevard Beaubourg 77184 Emerainville -

FRANCE

Pho : +33 (0)1 60 95 15 70 - Fax : +33 (0)1 60 77 41

[Sitemap](#) - [Legal notices](#)