

LKv 3910

MediLine



| | |
|-------------------------------------|-----------------------------|
| Gross/net capacity | 361/336 litre |
| Ext. dimensions in mm (w/d/h) | 597/615/1840 |
| Int. dimensions in mm (w/d/h) | 440/435/1635 |
| Energy consumption in 365 days | 309 kWh |
| Climate rating | 7 ¹ |
| Refrigerant | R 600a |
| Airborne acoustical noise emissions | 48 dB(A) |
| Connection Rating | 1.5 A |
| Voltage | 220-240 V~ |
| Cooling system | forced-air |
| Defrosting | automatic |
| Temperature range | +3°C to +16°C |
| Gradient * | 4.1°C |
| Max. fluctuation ** | 3.6°C |
| Housing material/colour | steel/white |
| Door/lid material | Steel |
| Liner Material | moulded polystyrol white |
| Type of Control | electronic controller |
| Temperature display | außen digital |
| Fault: warning signal | visual and audible |
| Interface | RS 485 |
| volt-free alarm contact | yes |
| Adjustable Shelving | 5 |
| Shelf material | grid shelves plastic-coated |
| Shelf loading | 45 kg |
| Handle | Ergonomic door handle |
| Lock | fitted |
| Self-closing door | yes |
| Door hinges | right, reversible |
| Insulation | 40 mm |
| Gross/net weight | 71 / 65 kg |
| EAN-Nr. | 9005382196854 |

1 Climate class / Climate Rätting: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C
3 +10°C to +25°C **4** +10°C to +30°C **5** +10°C to +40°C **7** +10°C to +35°C

* Gradient as defined in EN 60068-3: the difference between the warmest and coldest average measurements, increased by their expanded uncertainty, throughout the duration.

** Max. fluctuation as per EN 60068-3: greatest fluctuation value determined throughout the duration of measurements.