

CELL CULTURE DISHES




PRODUCT FEATURES

NEST cell culture dishes are ideal for all types of culturing where microscopic observation is required. Numeric marks on the bottom of each dish allow users to identify the location of cells.

- High clarity polystyrene
- Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment, excellent cell adherence
- Stackable for easy storage and handling
- Sterilized by E-Beam
- Non-pyrogenic



Easy-tear sterile package



Vents for gas exchange









Grip rings available for 100mm dishes



Stackable for easy storage and handling

Cell Culture Dishes

	Cat. No.	Type(mm)	Height	Cell Growth Area(cm ²)	Recommended Medium Volume(mL)	Sterile	/Pack	/Case
	706001	35	9.0	9	1.8-2.7	+	20	500
	705001	60	12.0	21	4.2-6.30	+	20	500
	704004	100	17.7	55	11-16.50	+	5	300
	704001	100	17.7	55	11-16.50	+	20	300
	704201	100	17.7	55	11-16.50	+	20	300
	715001	150	23.1	152	30.4-45.6	+	5	100