
CELL CULTURE FLASKS



Cell Culture Products

PRODUCT FEATURES

Total 16 different kinds of NEST cell culture flasks are available to meet your requirement, T25/T75/T175/T225, vent cap or plug seal cap, treated or Non-treated.

- Packaged in Zip bags
- High clarity, 100% virgin polystyrene
- Frosted writing and clear graduations
- Stackable
- Sterilized by E-Beam
- Vent caps (with breathable 0.22um membrane)
- Clear lot number for batch traceability
- Non-Pyrogenic, DNase/RNasefree




Vent caps with 0.22 um hydrophobic filters to eansure gas exchange without contamination



Yellow caps for non-treated flasks



Both vent caps and plug seal caps are available



T25/T75/T175/T225 are available

Tissue Culture Treated Cell Culture Flasks

Cat. No.	Cell Growth Area(cm²)	Volume	Recommende d Medium Volume(mL)	Сар Туре	Sterile	/Pack	/Case
707001	25	50	5-7.5	Plug Seal	+	10	200
708001	75	250	15-22.5	Plug Seal	+	5	100
709001	175	750	35-52.5	Plug Seal	+	5	40
721001	225	950	45-67.5	Plug Seal	+	5	25
707003	25	50	5-7.5	Vent	+	10	200
708003	75	250	15-22.5	Vent	+	5	100
709003	175	750	35-52.5	Vent	+	5	40
721003	225	950	45-67.5	Vent	+	5	25

Non-treated Cell Culture Flasks

	Cat. No.	Cell Growth Area(cm²)	Volume	Recommended Medium Volume(mL)	Сар Туре	Sterile	/Pack	/Case
MILL &	707011	25	50	5-7.5	Plug Seal	+	10	200
	708011	75	250	15-22.5	Plug Seal	+	5	100
	709011	175	750	35-52.5	Plug Seal	+	5	40
and the second s	721011	225	950	45-67.5	Plug Seal	+	5	25
Married Labor	707013	25	50	5-7.5	Vent	+	10	200
EIII	708013	75	250	15-22.5	Vent	+	5	100
(muutti	709013	175	750	35-52.5	Vent	+	5	40
The second second	721013	225	950	45-67.5	Vent	+	5	25
arrent Lake		1				1		