



CELL CULTURE PLATE



PRODUCT FEATURES

NEST cell culture plates have round edges for better handling. Low profile and low-evaporation lid design assures uniform growth. Both tissue culture treated and non-treated surfaces are available. Non-treated growth surface is more hydrophobic than tissue culture surface and is similar to the surface of a bacteriological Petri dish. Small feet on the bottom corners of each plate eliminate scraping of the microscope lens.


- Sterilized by E-beam
- High clarity, 100% virgin polystyrene
- Vacuum Plasma tissue culture treatment
- Clear lot number for batch traceability
- Peel-to-open, individually wrapped
- Markings of well coordinates available for 12, 24, 48, 96 and 384 well plates
- Non-Pyrogenic, DNase/RNase Free
- The package box of Nest cell culture plates can be used as reservoirs




Excellent bottom uniformity



Special lid design ensures low evaporation













Frosted and noise-free bottom design for easy grip












Lot No. for traceability

Tissue Culture Treated Plates

	Cat. No.	Type(well)	Cell Growth Area(cm ²)	Recommended Medium Volume(mL)	Colour	Bottom	Sterile	/Pack	/Case
	703001	6	9.5	1.90-2.90	Clear	Flat	+	1	50
	712001	12	3.6	0.76-1.14	Clear	Flat	+	1	50
	702001	24	1.9	0.38-0.57	Clear	Flat	+	1	50
	748001	48	0.88	0.19-0.285	Clear	Flat	+	1	50
	701001	96	0.32	0.10-0.20	Clear	Flat	+	1	100
	701002	96	0.32	0.10-0.20	Clear	Flat	+	10	100
	701101	96	-	0.10-0.20	Clear	U	+	1	100
	701201	96	-	0.10-0.20	Clear	V	+	1	100
	761001	384	0.11	0.025-0.05	Clear	Flat	+	1	100
	761601	384	0.11	0.025-0.05	White	Flat	+	1	100

Non-treated Cell Culture Plate

	Cat. No.	Type(well)	Max Well Volume(mL)	Colour	Bottom	Sterile	/Pack	/Case
	703011	6	16.8	Clear	Flat	+	1	50
	712011	12	6.9	Clear	Flat	+	1	50
	702011	24	3.4	Clear	Flat	+	1	50
	748011	48	1.6	Clear	Flat	+	1	50
	701011	96	0.36	Clear	Flat	+	1	100
	701111	96	0.33	Clear	U	+	1	100
	701211	96	0.32	Clear	V	+	1	100
	761011	384	0.11	Clear	Flat	+	1	100
	761611	384	0.11	White	Flat	+	1	100