

CELL SCRAPERS




Cell Culture Products

PRODUCT FEATURES

- Blade is made of HIPS and handle is made of PP
- Free rotating blade available to twist to requeired direction
- Two sizes are available
- Easy-tear package
- Non-pyrogenic and DNase/RNasefree
- Sterilized by E-beam, SAL=10⁻⁶







Cat. No.	Total Length(mm)	Blade Width(mm)	Material	Application	/Case
710001	220	13	HIPS/PP	T-25/T-75 Cell culture flask 60mm Cell culture dish	100
710011	280	20	HIPS/PP	T-75/T-175/T-225 Cell culture flask 100 and 150mm Cell culture dish	100