

CELL STRAINERS



Cell Culture Products

PRODUCT FEATURES

NEST cell strainer can be used to obtain uniform single cell suspension which is compatible with 50ml centrifuge tubes of most popular brands on the market.

- PP frame, Nylon strainer
- 3 sizes with different colors
- Sterile, individual package
- Non-pyrogenic and DNase/RNase free, non-cytotoxicity

Applicaton:

- Sample filtration
- Sample preparation for Flow Cytometer
- Isolating primary cells to consistently obtain a uniform single cell suspension from tissues

| Cat. No. | Filter diameter(μ L) | Mesh number | Color | Material | Sterile | /Pack | /Case |
|----------|---------------------------|-------------|--------|----------|---------|-------|-------|
| 258369 | 40 | 360 | Blue | PP/Nylon | + | 1 | 50 |
| 258368 | 70 | 220 | White | PP/Nylon | + | 1 | 50 |
| 258367 | 100 | 150 | Yellow | PP/Nylon | + | 1 | 50 |