

15/50 ML CENTRIFUGE TUBES



Sample Storage and Handling





PRODUCT FEATURES

- Printed clear graduations with marking area
- PP for the tube and HDPE for the cap
- High centrifugal strength up to 18,000 xg
- One-hand operation design
- Sterilized by E-Beam
- DNase, RNase and pyrogen free



Centrifuge Tubes, Premium

- Meet the requirements of analytical level
- High strength non-toxic materials
- Raw material comply with USP VI requirements

| | Cat. No. | Volume (mL) | Max.RCF(xg) | Packaging | Sterile | /Pack | /Case |
|---|----------|-------------|-------------|-----------|---------|-------|-------|
|  | 601001 | 15 | 15000 | PS Rack | + | 50 | 500 |
|  | 601002 | 15 | 15000 | Bulk | + | 50 | 500 |
|  | 602001 | 50 | 15000 | PS Rack | + | 25 | 500 |
|  | 602002 | 50 | 15000 | Bulk | + | 25 | 500 |



● Centrifuge Tubes, Standard


- Made from virgin polypropylene to meet the general requirements of the laboratory

● Centrifuge Tubes, Self-standing

- Made from virgin polypropylene, skirted and self-standing



| Cat. No. | Volume(mL) | Max.RCF(xg) | Packaging | Sterile | /Pack | /Case |
|----------|-------------------|-------------|-----------|---------|-------|-------|
| 601051 | 15 | 12000 | PS Rack | + | 50 | 500 |
| 601052 | 15 | 12000 | Bulk | + | 50 | 500 |
| 602051 | 50 | 12000 | PS Rack | + | 25 | 500 |
| 602052 | 50 | 12000 | Bulk | + | 25 | 500 |
| 602072 | 50(self-standing) | 5000 | Bulk | + | 25 | 500 |



● Centrifuge Tube Racks

- Made from virgin polypropylene, reusable and economic

| Cat. No. | Volume (mL) | Description | /Pack | /Case |
|----------|-------------|----------------------|-------|-------|
| 610001 | 15 | Individually wrapped | 1 | 50 |
| 610101 | 50 | Individually wrapped | 1 | 50 |