

## **Cover Glasses**



**Cell Culture Products** 

## **PRODUCT FEATURES**

- High quality glass, non-cytotoxic
- Thickness: 0.16-0.19mm
- · Manufactured in 100,000 grade clean room
- Resistance to organic solvent
- Autoclavable
- Can be used in confocal experiment

Cat. No.	Size(mm)	Application	Sterile	/Pack
801007	15	24 Well Cell Culture Plate	+	100
801008	20	12 Well Cell Culture Plate	+	100
801009	25	6 Well Cell Culture Plate	+	100

Tips

Coating

If the cover glass needs to be sterilized after long storage, please refer to the following autoclave condition: packed in tinfoil and put into a metal vessel,  $180^{\circ}\text{C}$ , dry-heating for 2 hours

It is suggested that the cover glass be coated with the following reagents to improve the cell adherence ability: polylysine, collagen, ECM or reagents with LN, fibronectin and polysaccharide