



# ELISA PLATES (IMMUNOPLATES)



Immunological Assay



# PRODUCT FEATURES

- PS material, designed for ELISA method
- Conform to SBS/ ANSI
- Both detachable and undetectable plates are available
- Flat and transparent bottom
- Excellent well uniformity, CV <5%



## ● Detachable ELISA Plate

- Clear alphanumeric codes for better test flexibility and identification
- Special designed stripe-holder to eliminate the change of disorientation
- 1\*8 stripes design for easy handling and readouts



Cat. No.	Description	/Pack	/Case
504201	96 Well ELISA Plate, Detachable, High Binding, Clear	5	50
514201	96 Well ELISA Plate, Undetectable, High Binding, Clear	5	50