




PCR® Tube Racks 96-Well

Perfect for pre- or post-PCR® sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes
- **Autoclavable**

Item No.	Color	Item No.	Color
HS23461A	Assorted* 	HS23461C	Blue 
HS23461B	Natural 	HS23461D	Green 

*Assorted colors (Blue, Green, Pink, Yellow, Orange)




Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
0.2	96	5.2	5 x 3.9 x 1.2	12.7 x 9.9 x 3.0	5 ea/pk, 40 pk/cs



PCR® Tube Racks Reversible

Removable hinged lid fits both sides

- Made of polypropylene
- PCR® side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Reverse side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes
- **Autoclavable**

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
HS2344A		0.2	168	5.5	8.1 x 5.2 x 2.1	20.6 x 13.1 x 5.3	5 ea/pk, 5 pk/cs
		0.5	40	7.6			
		1.5/2.0	12	11.4			

*Assorted colors (Blue, Green, Pink, Yellow, Orange)