

Rota-Rack® Tube Racks

Get organized with this fun, interlocked rack

- Made from polypropylene
- Modular design rotates and locks into place
- Alphanumeric identification for quick location and retrieval of tubes
- Suitable for a variety of culture, centrifuge, microcentrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tubes
- **Autoclavable**

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
HS29040A	Yellow	50	2	29	8.8 x 4.3 x 4.3	22.4 x 11 x 11	1 ea, 24 ea/cs
	Green	10-15	8	20			
	Blue	15-20	6	16			
	Purple	10-13	10	13.5			
HS29040B	Yellow	15	6	16.4	8.2 x 3.9 x 3.9	21 x 10 x 10	1 ea, 24 ea/cs
	Green	1.5/2.0	9	11.2			
	Blue	0.5	12	8			
	Purple	0.2	32	6.5			



Rotate and Customize
to suit your application needs

Rota-Rack® Duo Tube Rack

Dilution and concentration studies made easy

- Made from polypropylene
- Transfer solutions between 0.2 mL PCR® tubes and 0.5 and 1.5 mL microcentrifuge tubes
- Purple section holds 48 x 0.5 mL and 48 x 1.5/2.0 mL microfuge tubes on the reverse side
- Tapered wells securely suspend tubes and makes removal easy
- Blue section holds 96 x 0.2 mL tubes in a standard array pattern and 48 x 1.5/2.0 mL tubes on the reverse side
- Alphanumeric identification for quick location and retrieval of tubes
- **Autoclavable**

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm	UOM
HS29041A	Blue/ Purple	0.2	96	5	11.6 x 4.6 x 1.3	29.6 x 11.8 x 3.4	1 ea, 40 ea/cs
		0.5	48	8			
		1.5/2.0	96	11			

Easy Identification

Alphanumeric well identification

