

NEW!**Palm-shaped Lid**

suitable for left or right handed users
allows quick extraction of tubes

The lid can
be unlocked
and opened with
one hand.

Weighted Rubberized Base

eliminates potential movement and
further dampens vibration and noise

Sprout® & Sprout® Plus Mini Centrifuges

Ergonomic design | easy to operate | increased efficiency

With an ergonomic design and variety of colors, the Sprout® and Sprout® Plus Mini Centrifuges by Heathrow Scientific are loaded with features that bring increased efficiency to the lab.

Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use. The lid was functionally designed to fit in the palm of the hand for easy repeat opening and quick tube extraction. A weighted rubberized base eliminates potential movement and dampens vibration and noise.

Sprout® Plus

The Sprout Plus has the same footprint as the original Sprout, with a few added features and benefits.

- A locking lid for added protection when the unit is in use.
- Larger lid to accommodate longer and skirted tubes and taller tube lids
- Meets the latest applicable North American and EU safety requirements.

Sprout® Plus

Item No.	Color	L x W x H in	L x W x H cm	UOM
120610	Green	6 x 5 x 4.3	15.3 x 12.8 x 10.9	1 ea, 6 ea/cs
120611	Blue			

Sprout®

Item No.	Color	L x W x H in	L x W x H cm	UOM
HS120301	Green			
HS120395	Yellow	6 x 5 x 4.05	15.3 x 12.8 x 10.4	1 ea, 6 ea/cs
HS120396	Purple			

The Sprout and Sprout Plus are ideal for microfilter cell separations and HPLC samples, with adapters to accommodate all major microtube and PCR tube styles (1.5/2.0, 0.5, 0.2 mL) for customized use.

The Sprout and Sprout Plus centrifuge housings are made from ABS plastic and lids are polycarbonate plastic enabling the unit to easily cleaned.

A braking function quickly stops the unit when powered off or when the lid is opened.

**SnapSpin Rotors**

change out with the push of a button.
No tools are needed!

**Larger Lid**

on the Sprout Plus to accommodate skirted tubes and taller tube lids

**Specifications:**

Maximum Speed:	Approx. 6,000 rpm (fully loaded)
Maximum RCF:	2000 x g
Power:	100–240 VAC 50/60 HZ 0.6A, Approx. 7.3 W

Three-year warranty

**Package includes:**

- Sprout® or Sprout® Plus Mini-Centrifuge
- 6 place 1.5/2.0 mL tube rotor
- 16 place 0.2 mL tube rotor (singles, or 2 x 8 strip)
- 6 tube adapters 0.2 mL
- 6 tube adapters 0.5 mL
- 4 O-rings
- Storage case for rotors and adapters
- Tube rack
- Low voltage, double insulated power adapter with 4 interchangeable plugs





NEW!

Specially Designed Inserts
seat the centrifuge, all its accessories,
and includes a tube rack

12 Volt Car Adapter
for on-the-go testing

Portable Centrifuge Kit

Spin samples on-the-go

The Portable Centrifuge Kit by Heathrow Scientific offers a convenient way to spin down and prepare samples outside of a conventional laboratory environment.

The centrifuge kit enables the processing of samples before deterioration can take place and can be used in remote environments where a normal electrical supply is not an option. The kit enables immediate on-site testing for faster results.

The kit includes a Sprout® Mini Centrifuge, which comes with its own 12V car adapter and a 12V DC power adapter for when electrical outlets are in reach.

Conveniently carry the centrifuge and its components securely in the DuraPorter® transport box that has a specially designed insert to seat the centrifuge, all its accessories and includes a 16-place microtube rack. The DuraPorter is transparent, so any samples contained can be viewed and their integrity checked.

Package includes:

- Sprout® Mini Centrifuge
 - 6 Place 1.5/2.0 mL SnapSpin tube rotor
 - 16 Place 0.2 mL SnapSpin tube rotor (singles or 2x8 strip)
 - 6 Tube adapters (0.2 mL)
 - 6 Tube adapters (0.5 mL)
 - 4 O-rings
 - Storage case for rotors and adapters
 - 16 Place tube rack
 - Low voltage, double insulated power adapter with 4 interchangeable plugs



- DuraPorter® Transport Case
 - Tube rack insert (12 place 1.5/2.0 mL, 14 place 0.2 mL)
 - Transport box dividers
 - Sprout cradle
- Portable 12-volt Sprout Centrifuge Adapter

SEE PAGE 91 for more information on the DuraPorter®

Item No.		UOM
120568	Portable Centrifugation Kit, 100-240VAC, 50/60Hz Universal Plug + 12V Power Adapter	1 ea, 6 ea/cs
120569	Portable Centrifuge Transport Case with Inserts - only	1 ea
120570	12-Volt Sprout Power Adapter*	1 ea

*Adapter can be used with Sprout Mini Centrifuge units manufactured Jan 1, 2016 - present