

LGv 5010

MediLine



| | |
|-------------------------------------|-----------------------------|
| Gross/net capacity | 478/337 litre |
| Ext. dimensions in mm (w/d/h) | 747/750/1844 |
| Int. dimensions in mm (w/d/h) | 602/520/1224 |
| Energy consumption in 365 days | 1245 kWh |
| Climate rating | 7 ¹ |
| Refrigerant | R 290 |
| Airborne acoustical noise emissions | 55 dB(A) |
| Connection Rating | 3.0 A |
| Voltage | 220-240 V~ |
| Cooling system | forced-air |
| Defrosting | automatic |
| Temperature range | -9°C to -35°C |
| Gradient * | 5.8°C |
| Max. fluctuation ** | 4.9°C |
| Housing material/colour | steel/white |
| Door/lid material | Steel |
| Liner Material | moulded polystyrol white |
| Type of Control | electronic controller |
| Temperature display | außen digital |
| Fault: warning signal | visual and audible |
| Interface | RS 485 |
| volt-free alarm contact | yes |
| Number of drawers | 10 |
| Maximum number of baskets | 2 |
| Shelves | 5 |
| Shelf material | grid shelves plastic-coated |
| Shelf loading | 60 kg |
| Castors/type of castors | yes/swivel castors |
| Handle | Ergonomic door handle |
| Lock | fitted |
| Self-closing door | yes |
| Door hinges | right, reversible |
| Insulation | 70 mm |
| Defrost water drain | yes |
| Gross/net weight | 125 / 114 kg |
| EAN-Nr. | 9005382197394 |

1 Climate class / Climate Rätting: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C
3 +10°C to +25°C **4** +10°C to +30°C **5** +10°C to +40°C **7** +10°C to +35°C

* Gradient as defined in EN 60068-3: the difference between the warmest and coldest average measurements, increased by their expanded uncertainty, throughout the duration.

** Max. fluctuation as per EN 60068-3: greatest fluctuation value determined throughout the duration of measurements.