

8-STRIP PCR TUBES TUBES, CAPS



PRODUCT FEATURES

PCR Tube

- Made from USP grade VI polypropylene
- Compatible with most thermal cyclers and real-time PCR machine
- Extra-thin walls provide the optimal thermal transfer and shorter cycle times
- Easy open leak proof cap makes the loss of reaction less than 5%
- Free of Dnase and Rnase, human DNA and PCR inhibitor
- Autoclavable

| Cat. No. | Volume(mL) | Cap type | Color | /Pack | /Case |
|----------|------------|----------|-------|-------|--------|
| 401001 | 0.2 | Flat | Clear | 1,000 | 10,000 |



PCR 8-Strip Tubes

- NEST has both clear and white PCR strip tubes, which respectively apply to the general PCR and qPCR reactions
- Low profile tubes are designed to minimize the condensate contamination, particularly for the overnight experiment, and also reduce the evaporation of the reaction solution
- Free of Dnase and Rnase, human DNA and PCR inhibitor
- Autoclavable
- White strip tubes increase the fluorescent signal strength of real-time PCR



| Cat. No. | Volume(mL) | Specification | Color | /Pack | /Case |
|----------|------------|---------------|-------|-------|-------|
| 403102 | 0.1 | Low profile | Clear | 125 | 1250 |
| 403002 | 0.2 | High profile | Clear | 125 | 1250 |
| 403112 | 0.1 | Low profile | White | 125 | 1250 |
| 403012 | 0.2 | High profile | White | 125 | 1250 |



● PCR 8-Strip Tube with Flat Cap

- Eight individually capped thin wall tubes
- Caps are attached at an angle to prevent cap lids and hinges from interfering with each other
- Tight sealing, easy to open
- Flat, optically clear caps
- Extra-thin walls provide optimal thermal transfer and short cycle times
- Free of Dnase and Rnase, free of human DNA and PCR inhibitor
- Autoclavable

| Cat. No. | Volume(mL) | Specification | Color | /Pack | /Case |
|----------|------------|-------------------------------|-------|-------|-------|
| 404001 | 0.2 | High profile With flat cap | Clear | 120 | 1200 |



● PCR 8-Strip Tube Caps

- Domed and flat caps are available
- Compatible with both NEST PCR 8-strip tube and PCR plate
- Tight sealing, easy to open
- Flat cap is ideal for qPCR
- The design of domed cap can better prevent liquid evaporation
- Autoclavable

| Cat. No. | Specification | Color | /Pack | /Case |
|----------|---------------|-------|-------|-------|
| 406012 | Flat | Clear | 125 | 1250 |
| 406112 | Domed | Clear | 125 | 1250 |

