





PETRIDISHES




Bacterial Culture

PRODUCT FEATURES

- High clarity polystyrene
- Flat transparent surface
- Numeric mark on the bottom for easy identification
- PA/PE packaging
- Vent for gas exchange
- Stackable for easy storage and handling
- Sterilized by E-Beam, SAL10⁻⁶
- Non-pyrogenic





Contact dish

RODAC Dish is applied in hygiene monitoring to detect the possible source of contamination by contact cultures which is mostly used in pharmaceutical and food industries.

- Graduation into 10 mm squares with numbers to locate colonies
- Convex bottom designed to enable quick analysis
- Sterilized by E-Beam



Cat. No.	Description	Height(mm)	Volume(mm)	Sterile	/Pack	/Case
706011	35 mm	12	5	+	20	500
754001	60 mm	15	15	+	20	500
752001	90 mm	15	40	+	20	500
752011	90mm, 2 compartments	15	20*2	+	20	500
752021	90mm, 3 compartments	15	13*3	+	20	500
752031	90mm, 4 compartments	15	10*4	+	20	500
715011	150 mm	15	150	+	10	100
722011	60mm RODAC	15	15	+	20	500

