


RESERVOIR



Deep Well Plates and Reservoirs

PRODUCT FEATURES

- Low profile design, suitable for small volume robotic tips used in high flux instrument.
- In the automatic pipetting process, the reservoir can perfectly overcome the surface tension of the liquid and minimize residual liquid
- Matched reservoir breaker (for 96 and 384) and the cover can avoid the liquid wastage in the process of operation.
- Can be stored at -80°C
- Low heavy metal content
- Used high quality, imported PP material for high stability, ensuring no chemical reactions with test reagents.
- Compatible with DMSO and inert to water
- Certified DNase/RNase and Pyrogen Free
- Conformed to ANSI/SBS Standards



Cat. No.	Description	Volume(mL)	Working Volume(mL)	/Pack	/Case
360101	8-Channel, Reservoir, Without cap	22	-	5	50
360102	12-Channel, Reservoir, Without cap	12	-	5	50
360103	196-Channel, Reservoir, Without cap	195	-	5	50
360104	185-Channel, Reservoir, Without cap	185	-	5	50