SEROLOGICAL PIPETTES



Cell Culture Products

PRODUCT FEATURES

NEST serological pipettes can be used in both tissue culture and bacterial culture experiments.

- High clarity polystyrene
- Individually wrapped
- Edged color-coded band for identification
- Negative graduations for extravolume
- Sterilized by E-beam
- Non-pyrogenic, non-cytotoxic and non-haemolytic
- DNase/RNase free
- Polyolefin fiber filter to reduce contamination





Cat. No.	Size (mL)	Length(mm)	Color	/Box	/Case
324001	1	280	Yellow	500	3000
325001	2	280	Green	400	2400
326001	5	340	Blue	200	800
327001	10	339	Orange	200	800
328001	25	309	Red	200	800
329001	50	344	Purple	100	600