

gloveon

Alexa



Nitrile Exam Gloves Powder Free, Gauntlet Cuff

GloveOn® Alexa is the bold and powerful addition to our nitrile glove range, a tougher and stronger solution for workers in the healthcare sector. Distinguished by its gauntlet cuff that extends 500mm, GloveOn® Alexa is the glove of choice for maternity nurses assisting with water births, veterinarians, life science workers and room technicians responsible for cleaning surgical equipment.



		GloveOn® Alexa	
Length (mm)	500		
Thickness Measurements (mm)			
Palm (centre of Palm)	0.17 ± 0.02		
Finger (13mm ± 3mm from tip)	0.21 ± 0.02		
Physical Properties	Before Ageing	After Ageing	
Tensile Strength (MPa)	≥ 18	≥ 17	
Elongation (%)	≥ 500	≥ 400	
Inspection Levels & AQL	Inspection Level	AQL	
Watertightness	G1	1.5	
Physical Dimensions	S2	4.0	
Tensile Strength	S2	4.0	
Visual Inspection (Major)	S4	2.5	
Visual Inspection (Minor)	S4	4.0	
Particulate Residue	N = 5	≤ 2mg/glove	

REORDER CODE

NTR90SS SMALL
NTR90MM MEDIUM
NTR90LL LARGE
NTR90XL X-LARGE

FEATURES

- Palm textured
- Powder free
- Not made from natural rubber latex
- Ambidextrous
- Gauntlet cuff
- Maya blue colour

PACKAGING

50 gloves per box for S to L
45 gloves per box for XL
6 boxes per carton

REGULATORY COMPLIANCE

TGA – ARTG 280568, FDA 510(k),
MDD 93/42/EEC

STANDARDS

AS/NZS 4011.1:2014, EN 455 part 1, 2, 3,
HACCP Certified

MANUFACTURING ACCREDITATIONS

ISO 9001:2008
ISO 13485:2003

A brand by



Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com

