

gloveon Hammer

Nitrile Exam Gloves Powder Free, Standard Cuff



GloveOn® Hammer is the nitrile glove designed for the tough jobs that demand high tensile strength, elongation and a definitive colour. These gloves are suitable for security guards, tattoo artists, beauty stylists, motor mechanics, cleaners and general maintenance professionals.



GloveOn® Hammer		
Length (mm)	≥ 230	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.08 ± 0.02	
Finger (13mm ± 3mm from tip)	0.12 ± 0.02	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 500
Elongation after Isopropanol immersion (%)	≥ 500	N/A
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.5
Physical Dimensions	S2	4.0
Tensile Strength	S2	4.0
Visual Inspection (Major)	S4	2.5
Visual Inspection (Minor)	S4	4.0
Particulate Residue	N = 5	≤ 2mg/glove

REORDER CODE

NTR23XS	X-SMALL
NTR23SS	SMALL
NTR23MM	MEDIUM
NTR23LL	LARGE
NTR23XL	X-LARGE

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- UL certified to meet NFPA 1999
- Ambidextrous
- Black colour

PACKAGING

100 gloves per box for XS to L
90 gloves per box for XL
10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k),
MDD 93/42/EEC, REACH,
ROHS Directive 2002/95/EC

STANDARDS

ASTM D6319, ASTM D6978,
EN 455 part 1, 2, 3 & 4,
HACCP Certified, NFPA 1999

MANUFACTURING ACCREDITATIONS

ISO 9001:2015
ISO 13485:2003
EN ISO 13485:2012

A brand by

D:mun
www.munglobal.com

Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com

 **Birthing Kit
Foundation
Australia**
Proud Supporter of
Birthing Kit Foundation (Australia)

 Protection Always On