

gloveon

Matador



Nitrile Exam Gloves Powder Free, Long Cuff

GloveOn® Matador is built with quality, consistency and reliability in an 8ml nitrile glove. With GloveOn® Matador, users can take confidence in wearing a stronger glove to shape their performance in their work duties.



GloveOn® Matador		
Length (mm)	≥ 280	
Thickness Measurements (mm)		
Palm (centre of Palm)	0.20 ± 0.04	
Finger (13mm ± 3mm from tip)	0.28 ± 0.04	
Cuff (25mm from cuff edge)	0.14 ± 0.03	
Physical Properties	Before Ageing	After Ageing
Tensile Strength (MPa)	≥ 18	≥ 16
Elongation (%)	≥ 500	≥ 400
Inspection Levels & AQL	Inspection Level	AQL
Watertightness	G1	1.50
Physical Dimensions	S2	4.00
Tensile Strength	S2	4.00
Visual Inspection (Major)	S4	2.50
Visual Inspection (Minor)	S4	4.00
Particulate Residue	N = 5	≤ 2mg/glove

REORDER CODE

NTR18SS	SMALL
NTR18MM	MEDIUM
NTR18LL	LARGE
NTR18XL	X-LARGE
NTR18XXL	XX-LARGE

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Long cuff
- Sky blue colour

PACKAGING

50 gloves per box for S to XL
45 gloves per box for XXL
10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k), MDD 93/42/EEC, REACH, ROHS Directive 2002/95/EC, PPE 89/686/EEC

STANDARDS

ASTM D6319, EN 455 part 1, 2, 3 & 4, EN 374 part 1, 2, 3, AS/NZS 4011.1:2014

MANUFACTURING

ACCREDITATIONS

ISO 9001:2015
ISO 13485:2003
EN ISO 13485:2012

Tested in accordance with EN 374-3:2003 with the following chemical permeation performance levels:

Chemical	Performance Level	Permeation Performance Level	Measured Breakthrough Time (Minutes)
4% Chlorhexidine gluconate	6		
40% Sodium hydroxide	6		
35% Gluteraldehyde	6	1	> 10
1% Gluteraldehyde	6	2	> 30
4 % Gluteraldehyde	6	3	> 60
10% Formalin (Formaldehyde)	6	4	> 120
35% Ethanol	3	5	> 240
n-Hexane	2	6	> 480
Isopropanol	3		
Acetone	Not recommended		
3% Hydrogen Peroxide	6		
n-Heptane (J)	5		
96% Sulphuric Acid (K)	2		
40% Sodium Hydroxide (L)	6		

A brand by



Mun Australia
Toll free: 1800 456 837
munglobal.com.au



Proud Supporter of
Birthing Kit Foundation (Australia)



Protection Always On