

primeen

A brand by

mun

www.munglobal.com

en

Protection Always On



Shower Seal

Cast & Bandage Protector

Keep your cast dry while you shower

- Reusable
- Single patient use
- Latex free
- Made from heavy duty PVC
- Available in adult and child sizes

Shower Seal is the next step in protective wear for casts and bandages. The easy-to-use protection solution allows you to continue your normal showering routine, while keeping your plaster cast, dressing or bandage dry. Made from heavy duty PVC, Shower Seal is a reusable, reliable and quality waterproof protection solution for both adults and children.

Adult sizes				
Re-order code	CAA22SS-2	CAA39LL-2	CAL23SS-2	CAL42LL-2
Product	Adult Short Arm	Adult Long Arm	Adult Short Leg	Adult Long Leg
Seal Ring Diameter (mm)	150	200	200	250
Seal Shield Diameter (mm)	560	1000	580	1070
Carton QTY	20 Seals	20 Seals	20 Seals	20 Seals

Child sizes		
Re-order code	CCA18SS-2	CCL18SS-2
Product	Child Standard Arm	Child Standard Leg
Seal Ring Diameter (mm)	150	200
Seal Shield Diameter (mm)	460	460
Carton QTY	20 Seals	20 Seals

Fit Instructions

To determine the correct product length:

- Arm protector - measure from the tip of the finger to the top of the cast or bandage and add 2.5cm / 1".
- Leg protector - measure from the heel to the top of the cast or bandage and add 2.5cm / 1".

A brand by



Mun Australia
Toll free: 1800 456 837 (GLOVES)
info.au@munglobal.com



Protection Always On